Compte Rendu.

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



x10

x20

The sample consists mainly vital new bone (NB) with osteocytes and vascular structures. Areas of residual scaffold were surrounded by connective tissue and occupied ~10% of the sample area

NB- new bone
RS- residual scaffold
CT- connective tissue

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






RÉFÉRENCES

- (1) Ricci JL, Alexander H, Nadkarni P, Hawkins M, Turner J, Rosenblum SF, Brezenoff L, De Leonardis D, Pecora G. Biological mechanisms of calcium-sulfate replacement by bone. In: Davies JE, editor. *Bone Engineering*. Toronto: Em Squared Inc.; 2000. pp 332-344.
- (2) Anusavice KJ. Gypsum products. In: Anusavice KJ, editor. *Phillips' Science of Dental Materials*. St. Louis, MO: Saunders; 2003. pp 255-281.
- (3) Gao C, Huo S, Li X, You X, Zhang Y, Gao J. Characteristics of calcium sulfate/gelatin composite biomaterials for bone repair. *J Biomater Sci Polym Ed* 2007; 18: 799-824. [Links](#)
- (4) Neira-Carrillo A, Yazdani-Pedram M, Retuert J, Diaz-Dosque M, Gallois S, Arias JL. Selective crystallization of calcium salts by poly(acrylate)-grafted chitosan. *J Colloid Interface Sci* 2005; 286: 134-141. [Links](#)
- (5) Wright HB. Ridge Preservation Using Pre-Formed, Root-Shaped Calcium Sulfate Inserts. Lexington, KY: University of Kentucky College of Dentistry; 2004. 62p.
- (6) Guarnieri R, Grassi R, Ripari M, Pecora G. Maxillary sinus augmentation using granular calcium sulfate (surgiplaster sinus): Radiographic and histologic study at 2 years. *Int J Periodontics Restorative Dent* 2006; 26: 79-85. [Links](#)
- (7) Guarnieri R, Bovi M. Maxillary sinus augmentation using prehardened calcium sulfate: A case report. *Int J Periodontics Restorative Dent* 2002; 22: 503-508. [Links](#)
- (8) Orsini G, Ricci J, Scarano A, Pecora G, Petrone G, Iezzi G, Piattelli A. Bone-defect healing with calcium-sulfate particles and cement: An experimental study in rabbit. *J Biomed Mater Res B Appl Biomater* 2004; 68: 199-208. [Links](#)
- (9) Lioliou MG, Paraskeva CA, Koutsoukos PG, Payatakes AC. Calcium sulfate precipitation in the presence of watersoluble polymers. *J Colloid Interface Sci* 2006; 303: 164-170. [Links](#)
- (10) Wataha JC. Biocompatibility of dental materials. In: Anusavice KJ, editor. *Phillips' Science of Dental Materials*. St. Louis, MO: Saunders/Elsevier; 2003. pp 171-202.
- (11) Lillo R, Peltier LF. The substitution of plaster of Paris rods for portions of the diaphysis of the radius in dogs. *Surg Forum* 1956; 6: 556-558. [Links](#)
- (12) Alderman NE. Sterile plaster of Paris as an implant in the infrabony environment: A preliminary study. *J Periodontol* 1969; 40: 11-13. [Links](#)
- (13) Bahn SL. Plaster: A bone substitute. *Oral Surg Oral Med Oral Pathol* 1966; 21: 672-681. [Links](#)
- (14) Bell WH. Resorption characteristics of bone and bone substitutes. *Oral Surg Oral Med Oral Pathol* 1964; 17: 650-657. [Links](#)
- (15) Peltier LF, Bickel EY, Lillo R, Thein MS. The use of plaster of Paris to fill defects in bone. *Ann Surg* 1957; 146: 61-69. [Links](#)
- (16) Peltier LF. The use of plaster of Paris to fill defects in bone. *Clin Orthop* 1961; 21: 1-29. [Links](#)
- (17) Peltier LF, Orn D. The effect of the addition of plaster of Paris to autologous and homologous bone grafts in dogs. *Surg Forum* 1958; 8: 571-574. [Links](#)
- (18) Peltier LF. The use of plaster of Paris to fill large defects in bone: A preliminary report. *Am J Surg* 1959; 97: 311-315. [Links](#)
- (19) Shaffer CD, App GR. The use of plaster of Paris in treating infrabony periodontal defects in humans. *J Periodontol* 1971; 42: 685-690. [Links](#)
- (20) Kelly CM, Wilkins RM, Gitelis S, Hartjen C, Watson JT, Kim PT. The use of a surgical grade calcium sulfate as a bone graft substitute: Results of a multicenter trial. *Clin Orthop Relat Res* 2001; 42-50. [Links](#)
- (21) Scarano A, Orsini G, Pecora G, Iezzi G, Perrotti V, Piattelli A. Peri-implant bone regeneration with calcium sulfate: A light and transmission electron microscopy case report. *Implant Dent* 2007; 16: 195-203. [Links](#)
- (22) Wang ML, Massie J, Allen RT, Lee YP, Kim CW. Altered bioreactivity and limited osteoconductivity of calcium sulfate-based bone cements in the osteoporotic rat spine. *Spine J* 2007; 8: 349-350. [Links](#)
- (23) Wang ML, Massie J, Perry A, Garfin SR, Kim CW. A rat osteoporotic spine model for the evaluation of bioresorbable bone cements. *Spine J* 2007; 7: 466-474. [Links](#)
- (24) Payne JM, Cobb CM, Rapley JW, Killoy WJ, Spencer P. Migration of human gingival fibroblasts over guided tissue regeneration barrier materials. *J Periodontol* 1996; 67: 236-244. [Links](#)
- (25) Coetzee AS. Regeneration of bone in the presence of calcium sulfate. *Arch Otolaryngol* 1980; 106: 405-409. [Links](#)

- (26) Strocchi R, Orsini G, Iezzi G, Scarano A, Rubini C, Pecora G, Piattelli A. Bone regeneration with calcium sulfate: Evidence for increased angiogenesis in rabbits. *J Oral Implantol* 2002; 28: 273-278. [Links](#)
- (27) Walsh WR, Morberg P, Yu Y, Yang JL, Haggard W, Sheath PC, Svehla M, Bruce WJ. Response of a calcium sulfate bone graft substitute in a confined cancellous defect. *Clin Orthop* 2003: 228-236. [Links](#)
- (28) Park EK, Lee YE, Choi JY, Oh SH, Shin HI, Kim KH, Kim SY, Kim S. Cellular biocompatibility and stimulatory effects of calcium metaphosphate on osteoblastic differentiation of human bone marrow-derived stromal cells. *Biomaterials* 2004; 25: 3403-3411. [Links](#)
- (29) Mailland M, Waelchli R, Ruat M, Boddeke HG, Seuwen K. Stimulation of cell proliferation by calcium and a calcimimetic compound. *Endocrinology* 1997; 138: 3601-3605. [Links](#)
- (30) Yamauchi M, Yamaguchi T, Kaji H, Sugimoto T, Chihara K. Involvement of calcium-sensing receptor (CaR) in osteoblastic differentiation of mouse MC3T3-E1 cells. *Am J Physiol Endocrinol Metab* 2004; 3: E608-616. [Links](#)
- (31) Kameda T, Mano H, Yamada Y, Takai H, Amizuka N, Kobori M, Izumi N, Kawashima H, Ozawa H, Ikeda K, et al. Calcium-sensing receptor in mature osteoclasts, which are bone resorbing cells. *Biochem Biophys Res Commun* 1998; 245: 419-422. [Links](#)
- (32) Pi M, Hinson TK, Quarles L. Failure to detect the extracellular calcium-sensing receptor (CasR) in human osteoblast cell lines. *J Bone Miner Res* 1999; 14: 1310-1319. [Links](#)
- (33) Palmieri A, Pezzetti F, Brunelli G, Scapoli L, Muzio LL, Scarano A, Martinelli M, Carinci F. Calcium sulfate acts on the miRNA of MG63E osteoblast-like cells. *J Biomed Mater Res B Appl Biomater* 2008; 84: 369-374. [Links](#)
- (34) Lazary A, Balla B, Kosa JP, Bacsi K, Nagy Z, Takacs I, Varga PP, Speer G, Lakatos P. Effect of gypsum on proliferation and differentiation of MC3T3-E1 mouse osteoblastic cells. *Biomaterials* 2007; 28: 393-399. [Links](#)
- (35) Bier SJ, Sinensky MC. The versatility of calcium sulfate: Resolving periodontal challenges. *Compend Contin Educ Dent* 1999; 20: 655-61; Quiz 662. [Links](#)
- (36) Radentz WH, Collings CK. The implantation of plaster of paris in the alveolar process of the dog. *J Periodontol* 1965; 36: 357-364. [Links](#)
- (37) Sottosanti J. Calcium sulfate: A biodegradable and biocompatible barrier for guided tissue regeneration. *Compendium* 1992; 13: 226-228, 230, 232-234. [Links](#)
- (38) Sottosanti J. Calcium sulfate: An aid to periodontal, implant and restorative therapy. *J Calif Dent Assoc* 1992; 20: 45-46, 48, 52 Passim. [Links](#)
- (39) Selvig KA, Kersten BG, Wikesjo UM. Surgical treatment of intrabony periodontal defects using expanded polytetrafluoroethylene barrier membranes: Influence of defect configuration on healing response. *J Periodontol* 1993; 64: 730-733. [Links](#)
- (40) Wikesjo UM, Selvig KA. Periodontal wound healing and regeneration. *Periodontol* 2000 1999; 19: 21-39. [Links](#)
- (41) Sottosanti JS. Calcium sulfate is a safe, resorbable barrier adjunct to implant surgical procedures. *Dent Implantol Update* 1993; 4: 69-73. [Links](#)
- (42) Sottosanti JS. Aesthetic extractions with calcium sulfate and the principles of guided tissue regeneration. *Pract Periodontics Aesthet Dent* 1993; 5: 61-69; Quiz 69. [Links](#)
- (43) Sottosanti JS. Calcium sulfate-aided bone regeneration: A case report. *Periodontal Clin Investig* 1995; 17: 10-15. [Links](#)
- (44) Sottosanti JS. Calcium sulfate: A valuable addition to the implant/bone regeneration complex. *Dent Implantol Update* 1997; 8: 25-29. [Links](#)
- (45) Kim CK, Kim HY, Chai JK, Cho KS, Moon IS, Choi SH, Sottosanti JS, Wikesjo UM. Effect of a calcium sulfate implant with calcium sulfate barrier on periodontal healing in 3-wall intrabony defects in dogs. *J Periodontol* 1998; 69: 982-988. [Links](#)
- (46) Kim CK, Chai JK, Cho KS, Moon IS, Choi SH, Sottosanti JS, Wikesjo UM. Periodontal repair in intrabony defects treated with a calcium sulfate implant and calcium sulfate barrier. *J Periodontol* 1998; 69: 1317-1324. [Links](#)
- (47) Rosen PS, Reynolds MA. Polymer-assisted regenerative therapy: Case reports of 22 consecutively treated periodontal defects with a novel combined surgical approach. *J Periodontol* 1999; 70: 554-561. [Links](#)
- (48) MacNeill SR, Cobb CM, Rapley JW, Glaros AG, Spencer P. In vivo comparison of synthetic osseous graft materials. A preliminary study. *J Clin Periodontol* 1999; 26: 239-245. [Links](#)
- (49) Cui X, Zhang B, Wang Y, Gao Y. Effects of chitosan-coated pressed calcium sulfate pellet combined with recombinant human bone morphogenetic protein 2 on restoration of segmental bone defect. *J Craniofac Surg* 2008; 19: 459-465. [Links](#)

- (⁵⁰) Melo LG, Nagata MJ, Bosco AF, Ribeiro LL, Leite CM. Bone healing in surgically created defects treated with either bioactive glass particles, a calcium sulfate barrier, or a combination of both materials. A histological and histometric study in rat tibias. *Clin Oral Implants Res* 2005; 16: 683-691. [Links](#)
- (⁵¹) Lewandrowski KU, Cattaneo MV, Gresser JD, Wise DL, White RL, Bonassar L, Trantolo DJ. Effect of a poly(propylenefumarate) foaming cement on the healing of bone defects. *Tissue Eng* 1999; 5: 305-316. [Links](#)
- (⁵²) Hing KA, Wilson LF, Buckland T. Comparative performance of three ceramic bone graft substitutes. *Spine J* 2007; 7: 475-490. [Links](#)
- (⁵³) Anson D. Calcium sulfate: A 4-year observation of its use as a resorbable barrier in guided tissue regeneration of periodontal defects. *Compend Contin Educ Dent* 1996; 17: 895-899. [Links](#)
- (⁵⁴) Anson D. Using calcium sulfate in guided tissue regeneration: a recipe for success. *Compend Contin Educ Dent* 2000; 21: 365-370, 372-373, 376; Quiz 378. [Links](#)
- (⁵⁵) Deliberador TM, Nagata MJ, Furlaneto FA, Melo LG, Okamoto T, Sundefeld ML, Fucini SE. Autogenous bone graft with or without a calcium sulfate barrier in the treatment of Class II furcation defects: A histologic and histometric study in dogs. *J Periodontol* 2006; 77: 780-789. [Links](#)
- (⁵⁶) Aichelmann-Reidy ME, Heath CD, Reynolds MA. Clinical evaluation of calcium sulfate in combination with demineralized freeze-dried bone allograft for the treatment of human intraosseous defects. *J Periodontol* 2004; 75: 340-347. [Links](#)
- (⁵⁷) Couri CJ, Maze GI, Hinkson DW, Collins BH III, Dawson DV. Medical grade calcium sulfate hemihydrate versus expanded polytetrafluoroethylene in the treatment of mandibular class II furcations. *J Periodontol* 2002; 73: 1352-1359. [Links](#)
- (⁵⁸) Orsini M, Orsini G, Benlloch D, Aranda JJ, Lazaro P, Sanz M, De Luca M, Piattelli A. Comparison of calcium sulfate and autogenous bone graft to bioabsorbable membranes plus autogenous bone graft in the treatment of intrabony periodontal defects: A split-mouth study. *J Periodontol* 2001; 72: 296-302. [Links](#)
- (⁵⁹) Setya AB, Bissada NF. Clinical evaluation of the use of calcium sulfate in regenerative periodontal surgery for the treatment of Class III furcation involvement. *Periodontal Clin Investig* 1999; 21: 5-14. [Links](#)
- (⁶⁰) Pecora G, Andreana S, Margarone JE III, Covani U, Sottosanti JS. Bone regeneration with a calcium sulfate barrier. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 1997; 84: 424-429. [Links](#)
- (⁶¹) Pecora GE, De Leonardis D, Della Rocca C, Cornelini R, Cortesini C. Short-term healing following the use of calcium sulfate as a grafting material for sinus augmentation: A clinical report. *Int J Oral Maxillofac Implants* 1998; 13: 866-873. [Links](#)
- (⁶²) De Leonardis D, Pecora GE. Augmentation of the maxillary sinus with calcium sulfate: One-year clinical report from a prospective longitudinal study. *Int J Oral Maxillofac Implants* 1999; 14: 869-878. [Links](#)
- (⁶³) De Leonardis D, Pecora GE. Prospective study on the augmentation of the maxillary sinus with calcium sulfate: Histological results. *J Periodontol* 2000; 71: 940-947. [Links](#)
- (⁶⁴) Tadjoeidin ES, de Lange GL, Holzmann PJ, Kulper L, Burger EH. Histological observations on biopsies harvested following sinus floor elevation using a bioactive glass material of narrow size range. *Clin Oral Implants Res* 2000; 11: 334-344. [Links](#)
- (⁶⁵) Tadjoeidin ES, de Lange GL, Bronckers AL, Lyaruu DM, Burger EH. Deproteinized cancellous bovine bone (Bio-Oss) as bone substitute for sinus floor elevation. A retrospective, histomorphometrical study of five cases. *J Clin Periodontol* 2003; 30: 261-270. [Links](#)
- (⁶⁶) Hallman M, Sennerby L, Lundgren S. A clinical and histologic evaluation of implant integration in the posterior maxilla after sinus floor augmentation with autogenous bone, bovine hydroxyapatite, or a 20:80 mixture. *Int J Oral Maxillofac Implants* 2002; 17: 635-643. [Links](#)
- (⁶⁷) Wetzel AC, Stich H, Caffesse RG. Bone apposition onto oral implants in the sinus area filled with different grafting materials. A histological study in beagle dogs. *Clin Oral Implants Res* 1995; 6: 155-163. [Links](#)
- (⁶⁸) Iezzi G, Fiera E, Scarano A, Pecora G, Piattelli A. Histologic evaluation of a provisional implant retrieved from man 7 months after placement in a sinus augmented with calcium sulphate: A case report. *J Oral Implantol* 2007; 33: 89-95. [Links](#)

- (⁶⁹) Scarano A, Degidi M, Iezzi G, Pecora G, Piattelli M, Orsini G, Caputi S, Perrotti V, Mangano C, Piattelli A. Maxillary sinus augmentation with different biomaterials: A comparative histologic and histomorphometric study in man. *Implant Dent* 2006; 15: 197-207. [Links](#)
- (⁷⁰) Andreana S, Cornelini R, Edsberg LE, Natiella JR. Maxillary sinus elevation for implant placement using calcium sulfate with and without DFDBA: Six cases. *Implant Dent* 2004; 13: 270-277. [Links](#)
- (⁷¹) Shi B, Zhou Y, Wang YN, Cheng XR. Alveolar ridge preservation prior to implant placement with surgical-grade calcium sulfate and platelet-rich plasma: A pilot study in a canine model. *Int J Oral Maxillofac Implants* 2007; 22: 656-665. [Links](#)
- (⁷²) Rolph CW. Ridge Preservation Using Bioactive Glass and a Calcium Sulfate Barrier. Lexington, KY: University of Kentucky College of Dentistry; 2000. 62p.
- (⁷³) Camargo PM, Lekovic V, Weinlaender M, Klokkevold PR, Kenney EB, Dimitrijevic B, Nedic M, Jancovic S, Orsini M. Influence of bioactive glass on changes in alveolar process dimensions after exodontia. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2000; 90: 581-586. [Links](#)
- (⁷⁴) Vance GS, Greenwell H, Miller RL, Hill M, Johnston H, Scheetz JP. Comparison of an allograft in an experimental putty carrier and a bovine-derived xenograft used in ridge preservation: A clinical and histologic study in humans. *Int J Oral Maxillofac Implants* 2004; 19: 491-497. [Links](#)
- (⁷⁵) Rosenblum SF, Frenkel S, Ricci JR, Alexander H. Diffusion of fibroblast growth factor from a plaster of Paris carrier. *J Appl Biomater* 1993; 4: 67-72. [Links](#)
- (⁷⁶) Nyan M, Sato D, Oda M, Machida T, Kobayashi H, Nakamura T, Kasugai S. Bone formation with the combination of simvastatin and calcium sulfate in critical-sized rat calvarial defect. *J Pharmacol Sci* 2007; 104: 384-386. [Links](#)
- (⁷⁷) Intini GE III. Engineering a bioactive matrix by modifications of calcium sulfate. *Tissue Eng* 2002; 8: 997-1008. [Links](#)
- (⁷⁸) Intini G, Andreana S, Intini FE, Buhite RJ, Bobek LA. Calcium sulfate and platelet-rich plasma make a novel osteoinductive biomaterial for bone regeneration. *J Transl Med* 2007; 5: 13. [Links](#)
- (⁷⁹) Bateman J, Intini G, Margarone J, Goodloe S, Bush P, Lynch SE, Dziak R. Platelet-derived growth factor enhancement of two alloplastic bone matrices. *J Periodontol* 2005; 76: 1833-1841. [Links](#)
- (⁸⁰) Kim SG, Chung CH, Kim YK, Park JC, Lim SC. Use of particulate dentin-plaster of Paris combination with/without platelet-rich plasma in the treatment of bone defects around implants. *Int J Oral Maxillofac Implants* 2002; 17: 86-94. [Links](#)
- (⁸¹) Frenkel SR, Simon J, Alexander H, Dennis M, Ricci JL. Osseointegration on metallic implant surfaces: Effects of microgeometry and growth factor treatment. *J Biomed Mater Res* 2002; 63: 706-713. [Links](#)
- (⁸²) Damien CJ, Parsons JR, Benedict JJ, Weisman DS. Investigation of a hydroxyapatite and calcium sulfate composite supplemented with an osteoinductive factor. *J Biomed Mater Res* 1990; 24: 639-654. [Links](#)
- (⁸³) Kennedy JE. Local Delivery of Antibiotics in a Calcium Sulfate Vehicle: In Vitro Results. Lexington, KY: University of Kentucky College of Dentistry; 2000. 63p.
- (⁸⁴) Doadrio JC, Arcos D, Cabanas MV, Vallet-Regi M. Calcium sulphate-based cements containing cephalexin. *Biomaterials* 2004; 25: 2629-2635. [Links](#)
- (⁸⁵) Richelsoph KC, Webb ND, Haggard WO. Elution behavior of daptomycin-loaded calcium sulfate pellets: A preliminary study. *Clin Orthop Relat Res* 2007; 461: 68-73. [Links](#)
- (⁸⁶) Webb ND, McCanless JD, Courtney HS, Bumgardner JD, Haggard WO. Daptomycin eluted from calcium sulfate appears effective against *Staphylococcus*. *Clin Orthop Relat Res* 2008; 466: 1383-1387. [Links](#)
- (⁸⁷) Yarboro SR, Baum EJ, Dahners LE. Locally administered antibiotics for prophylaxis against surgical wound infection. An in vivo study. *J Bone Joint Surg Am* 2007; 89: 929-933. [Links](#)
- (⁸⁸) Pepelassi EM, Bissada NF, Greenwell H, Farah CF. Doxycycline-tricalcium phosphate composite graft facilitates osseous healing in advanced periodontal furcation defects. *J Periodontol* 1991; 62: 106-115. [Links](#)
- (⁸⁹) DiBattista P, Bissada NF, Ricchetti PA. Comparative effectiveness of various regenerative modalities for the treatment of localized juvenile periodontitis. *J Periodontol* 1995; 66: 673-678. [Links](#)

- ⁽⁹⁰⁾ Mousset B, Benoit MA, Delloye C, Bouillet R, Gillard J. Biodegradable implants for potential use in bone infection. An in vitro study of antibiotic-loaded calcium sulphate. *Int Orthop* 1995; 19: 157-161. [Links](#)
- ⁽⁹¹⁾ Benoit MA, Mousset B, Delloye C, Bouillet R, Gillard J. Antibiotic-loaded plaster of Paris implants coated with polylactide-co-glycolide as a controlled release delivery system for the treatment of bone infections. *Int Orthop* 1997; 21: 403-408. [Links](#)
- ⁽⁹²⁾ Bibbo C, Patel DV. The effect of demineralized bone matrix-calcium sulfate with vancomycin on calcaneal fracture healing and infection rates: A prospective study. *Foot Ankle Int* 2006; 27: 487-493. [Links](#)
- ⁽⁹³⁾ Chang W, Colangeli M, Colangeli S, Di Bella C, Gozzi E, Donati D. Adult osteomyelitis: Debridement versus debridement plus Osteoset T pellets. *Acta Orthop Belg* 2007; 73: 238-243. [Links](#)
- ⁽⁹⁴⁾ Salgami EV, Bowling FL, Whitehouse RW, Boulton AJ. Use of tobramycin-impregnated calcium sulphate pellets in addition to oral antibiotics: An alternative treatment to minor amputation in a case of diabetic foot osteomyelitis. *Diabetes Care* 2007; 30: 181-182. [Links](#)
- ⁽⁹⁵⁾ Tuzuner T, Uygun I, Sencan I, Haklar U, Oktas B, Ozdemir D. Elution characteristics and mechanical properties of calcium sulfate-loaded bone cement containing teicoplanin. *J Orthop Sci* 2007; 12: 170-177. [Links](#)
- ⁽⁹⁶⁾ Tuzuner T, Sencan I, Ozdemir D, Alper M, Duman S, Yavuz T, Yildirim M. In vivo evaluation of teicoplanin- and calcium sulfate-loaded PMMA bone cement in preventing implant-related osteomyelitis in rats. *J Chemother* 2006; 18: 628-633. [Links](#)
- ⁽⁹⁷⁾ Vechasilp J, Tangtrakulwanich B, Oungbho K, Yuenyongsawad S. The efficacy of methotrexate-impregnated hydroxyapatite composites on human mammary carcinoma cells. *J Orthop Surg (Hong Kong)* 2007; 15: 56-61. [Links](#)
- ⁽⁹⁸⁾ Buranapanitkit B, Srinilta V, Ingviga N, Oungbho K, Geater A, Ovatlarnporn C. The efficacy of a hydroxyapatite composite as a biodegradable antibiotic delivery system. *Clin Orthop Relat Res* 2004: 244-252. [Links](#)
- ⁽⁹⁹⁾ Kim SY, Kim SG, Lim SC, Bae CS. Effects on bone formation in ovariectomized rats after implantation of tooth ash and plaster of Paris mixture. *J Oral Maxillofac Surg* 2004; 62: 852-857. [Links](#)
- ⁽¹⁰⁰⁾ Parsons JR, Ricci JL, Alexander H, Bajpai PK. Osteoconductive composite grouts for orthopedic use. *Ann N Y Acad Sci* 1988; 523: 190-207. [Links](#)
- ⁽¹⁰¹⁾ Ziran BH, Smith WR, Morgan SJ. Use of calcium-based demineralized bone matrix/allograft for nonunions and posttraumatic reconstruction of the appendicular skeleton: Preliminary results and complications. *J Trauma* 2007; 63: 1324-1328. [Links](#)
- ⁽¹⁰²⁾ Robinson D, Alk D, Sandbank J, Farber R, Halperin N. Inflammatory reactions associated with a calcium sulfate bone substitute. *Ann Transplant* 1999; 4(3-4): 91-97. [Links](#)
- ⁽¹⁰³⁾ Borrelli J Jr, Prickett WD, Ricci WM. Treatment of nonunions and osseous defects with bone graft and calcium sulfate. *Clin Orthop Relat Res* 2003; 411: 245-254. [Links](#)
- ⁽¹⁰⁴⁾ Furlaneto FA, Nagata MJ, Fucini SE, Deliberador TM, Okamoto T, Messori MR. Bone healing in critical-size defects treated with bioactive glass/calcium sulfate: A histologic and histometric study in rat calvaria. *Clin Oral Implants Res* 2007; 18: 311-318. [Links](#)
- ⁽¹⁰⁵⁾ Lewis KN, Thomas MV, Puleo DA. Mechanical and degradation behavior of polymer-calcium sulfate composites. *J Mater Sci Mater Med* 2006; 17: 531-537. [Links](#)
- ⁽¹⁰⁶⁾ Urban RM, Turner TM, Hall DJ, Infanger SI, Cheema N, Lim TH, Richelsoph K. An injectable calcium sulfate-based bone graft putty using hydroxypropylmethylcellulose as the plasticizer. *Orthopedics* 2004; 27(1 Suppl): s155-s159. [Links](#)
- ⁽¹⁰⁷⁾ Reynolds MA, Aichelmann-Reidy ME, Kassolis JD, Prasad HS, Rohrer MD. Calcium sulfate-carboxymethylcellulose bone graft binder: Histologic and morphometric evaluation in a critical size defect. *J Biomed Mater Res B Appl. Biomater* 2007; 83: 451-458. [Links](#)
- ⁽¹⁰⁸⁾ d'Ayala GG, De Rosa A, Laurienzo P, Malinconico M. Development of a new calcium sulphate-based composite using alginate and chemically modified chitosan for bone regeneration. *J Biomed Mater Res A* 2007; 81: 811-820. [Links](#)
- ⁽¹⁰⁹⁾ Krisanapiboon A, Buranapanitkit B, Oungbho K. Biocompatibility of hydroxyapatite composite as a local drug delivery system. *J Orthop Surg (Hong Kong)* 2006; 14: 315-318. [Links](#)
- ⁽¹¹⁰⁾ Mamidwar SS, Arena C, Kelly S, Alexander H, Ricci J. In vitro characterization of a calcium sulfate/PLLA composite for use as a bone graft material. *J Biomed Mater Res B Appl Biomater* 2007; 81: 57-65. [Links](#)

- (¹¹¹) Mamidwar S, Weiner M, Alexander H, Ricci J. In vivo bone response to calcium sulfate/poly L-lactic acid composite. *Implant Dent* 2008; 17: 208-216. [Links](#)
- (¹¹²) Huan Z, Chang J. Self-setting properties and in vitro bioactivity of calcium sulfate hemihydrate-tricalcium silicate composite bone cements. *Acta Biomater* 2007; 3: 952-960. [Links](#)
- (¹¹³) Urban RM, Turner TM, Hall DJ, Inoue N, Gitelis S. Increased bone formation using calcium sulfate-calcium phosphate composite graft. *Clin Orthop Relat Res* 2007; 459: 110-117. [Links](#)
- (¹¹⁴) Guo H, Wei J, Liu CS. Development of a degradable cement of calcium phosphate and calcium sulfate composite for bone reconstruction. *Biomed Mater* 2006; 1: 193-197. [Links](#)
- (¹¹⁵) Kim JY, Cho BC. Effect of calcium sulfate pellets on early bone mineralization in distraction osteogenesis for craniofacial microsomnia in adults. *J Craniofac Surg* 2007; 18: 1309-1321. [Links](#)
- (¹¹⁶) Cho BC, Park JW, Baik BS, Kim IS. Clinical application of injectable calcium sulfate on early bony consolidation in distraction osteogenesis for the treatment of craniofacial microsomnia. *J Craniofac Surg* 2002; 13: 465-475; Discussion 475-477. [Links](#)
- (¹¹⁷) Cho BC, Chung HY, Lee DG, Yang JD, Park JW, Roh KH, Kim GU, Lee DS, Kwon IC, Bae EH, et al. The effect of chitosan bead encapsulating calcium sulfate as an injectable bone substitute on consolidation in the mandibular distraction osteogenesis of a dog model. *J Oral Maxillofac Surg* 2005; 63: 1753-1764. [Links](#)
- (¹¹⁸) Cho BC, Kim TG, Yang JD, Chung HY, Park JW, Kwon IC, Roh KH, Chung HS, Lee DS, Park NU, Kim IS. Effect of calcium sulfate-chitosan composite: Pellet on bone formation in bone defect. *J Craniofac Surg* 2005; 16: 213-224; Discussion 225-227. [Links](#)
- (¹¹⁹) Song HR, Oh CW, Kyung HS, Park IH, Kim PT, Baek SH, Kim SJ, Lee ST. Injected calcium sulfate for consolidation of distraction osteogenesis in rabbit tibia. *J Pediatr Orthop B* 2004; 13: 170-175. [Links](#)
- (¹²⁰) Cho BC, Park JW, Baik BS, Kwon IC, Kim IS. The role of hyaluronic acid, chitosan, and calcium sulfate and their combined effect on early bony consolidation in distraction osteogenesis of a canine model. *J Craniofac Surg* 2002; 13: 783-793. [Links](#)
- (¹²¹) Mittal M, Chandra S, Chandra S. An evaluation of plaster of Paris barriers used under various materials to repair furcation perforations (in vitro study). *J Endod* 1999; 25: 385-388. [Links](#)
- (¹²²) Alhadainy HA, Abdalla AI. Artificial floor technique used for the repair of furcation perforations: A microleakage study. *J Endod* 1998; 24: 33-35. [Links](#)
- (¹²³) Conner HD. Bone grafting with a calcium sulfate barrier after root amputation. *Compend Contin Educ Dent* 1996; 17: 42, 44, 46; Quiz 48. [Links](#)
- (¹²⁴) Andreana S. A combined approach for treatment of developmental groove associated periodontal defect. A case report. *J Periodontol* 1998; 69: 601-607. [Links](#)
- (¹²⁵) Pecora G, De Leonardis D, Ibrahim N, Bovi M, Cornelini R. The use of calcium sulphate in the surgical treatment of a through and through periradicular lesion. *Int Endod J* 2001; 34: 189-197. [Links](#)
- (¹²⁶) Pecora G, Baek SH, Rethnam S, Kim S. Barrier membrane techniques in endodontic microsurgery. *Dent Clin North Am* 1997; 41: 585-602. [Links](#)
- (¹²⁷) Doobrow JH, Leite RS, Hirsch HZ. Concomitant oroantral communication repair and immediate implant placement: A five-year case report. *Implant Dent* 2008; 17: 176-181. [Links](#)
- (¹²⁸) Mesimeris V, Sade E, Baer PN. Calcium sulfate as a biodegradable barrier membrane - A preliminary report on the Surgiplast. Technique. *Periodontal Clin Investig* 1995; 17: 13-16. [Links](#)
- (¹²⁹) Anson D. Saving periodontally hopeless teeth using calcium sulfate and demineralized freeze-dried bone allograft. *Compend Contin Educ Dent* 1998; 19: 284, 286, 288 *Passim*. [Links](#)
- (¹³⁰) Chen WJ, Tsai TT, Chen LH, Niu CC, Lai PL, Fu TS, McCarthy K. The fusion rate of calcium sulfate with local autograft bone compared with autologous iliac bone graft for instrumented short-segment spinal fusion. *Spine* 2005; 30: 2293-2297. [Links](#)
- (¹³¹) Alexander DI, Manson NA, Mitchell MJ. Efficacy of calcium sulfate plus decompression bone in lumbar and lumbosacral spinal fusion: Preliminary results in 40 patients. *Can J Surg* 2001; 44: 262-266. [Links](#)
- (¹³²) Peters CL, Hines JL, Bachus KN, Craig MA, Bloebaum RD. Biological effects of calcium sulfate as a bone graft substitute in ovine metaphyseal defects. *J Biomed Mater Res A* 2006; 76: 456-462. [Links](#)

- (¹³³) Anson D. Calcium sulfate-augmented soft tissue root coverage adjacent to connective tissue grafting: A new technique. *Int J Periodontics Restorative Dent* 2003; 23: 337-343. Links
- (¹³⁴) Frame JW. A composite of porous calcium sulphate dihydrate and cyanoacrylate as a substitute for autogenous bone. *J Oral Surg* 1980; 38: 251-256. Links
- (¹³⁵) Perry A, Mahar A, Massie J, Arrieta N, Garfin S, Kim C. Biomechanical evaluation of kyphoplasty with calcium sulfate cement in a cadaveric osteoporotic vertebral compression fracture model. *Spine J* 2005; 5: 489-493. Links
- (¹³⁶) Derincek A, Wu C, Mehbod A, Transfeldt EE. Biomechanical comparison of anatomic trajectory pedicle screw versus injectable calcium sulfate graft-augmented pedicle screw for salvage in cadaveric thoracic bone. *J Spinal Disord Tech* 2006; 19: 286-291. Links
- (¹³⁷) Hak DJ. The use of osteoconductive bone graft substitutes in orthopaedic trauma. *J Am Acad Orthop Surg* 2007; 15: 525-536. Links
- (¹³⁸) Panchbhavi VK, Vallurupalli S, Morris R, Patterson R. The use of calcium sulfate and calcium phosphate composite graft to augment screw purchase in osteoporotic ankles. *Foot Ankle Int* 2008; 29: 593-600. Links
- (¹³⁹) Kelly CM, Wilkins RM. Treatment of benign bone lesions with an injectable calcium sulfate-based bone graft substitute. *Orthopedics* 2004; 27(1 Suppl): s131-s135. Links
- (¹⁴⁰) Wikesjo UM, Susin C, Qahash M, Polimeni G, Leknes KN, Shanaman RH, Prasad HS, Rohrer MD, Hall J. The critical size supraalveolar peri-implant defect model: characteristics and use. *J Clin Periodontol* 2006; 33: 846-854. Links
- (¹⁴¹) Raiche AT, Puleo DA. In vitro effects of combined and sequential delivery of two bone growth factors. *Biomaterials* 2004; 25: 677-685. Links
- (¹⁴²) Raiche AT, Puleo DA. Modulated release of bioactive protein from multilayered blended PLGA coatings. *Int J Pharm* 2006; 311(1-2): 40-49. Links
- (¹⁴³