

Calcium sulfate acts on stem cells derived from peripheral blood

Vincenzo Sollazzo,¹ Annalisa Palmieri,² Luca Scapoli,³
 Marcella Martinelli,³ Ambra Girardi,³ Furio Pezzetti,³
 Francesca Farinella,² Francesco Carinci^{2*}

Calcium sulfate (CaS) is a highly biocompatible material and enhances bone formation in vivo. However, how CaS alters osteoblast activity to promote bone formation is incompletely understood.

To study how CaS can induce osteoblast differentiation in mesenchymal stem cells, the expression levels of bone related genes and mesenchymal stem cell markers were analyzed, using real time Reverse Transcription-Polymerase Chain Reaction.

CaS causes a significant induction of the bone related genes osteopontin (SPP1), osteocalcin (BGLAP) and collagen, type III, alpha 1 (COL3A1). The expression of ENG and FOSL1 were not significantly changed in stem cells treated with CaS respect to untreated cells, while RUNX1, COL1A1 and ALPL were significantly down expressed.

The results obtained lend to a better understand of the molecular mechanism of bone regeneration and as a model for comparing other materials with similar clinical effects.

(J Osteol Biomat 2010; 1:55-60)

Key words: calcium sulfate, gene expression, stem cells, bone regeneration.

¹Orthopedic Clinic, University of Ferrara, Ferrara, Italy;

²Department of Maxillofacial Surgery, University of Ferrara, Ferrara, Italy;

³Department of Histology, Embryology and Applied Biology, University of Bologna, Bologna, Italy;”

Corresponding author:

*Prof. Francesco Carinci, MD

Department of D.M.C.C.C.

University of Ferrara

Corso Giovecca, 203 - 44100 Ferrara (Italy)

Phone/Fax : 0039-0532-455582

E-mail: crc@unife.it - Web: www.carinci.org

INTRODUCTION

Several graft materials have been proposed in implants dentistry. Autogenous bone is the gold standard but usually donor oral sites have a limited amount of graft material. Consequently, surgeons harvest bone from extra-oral sites with increased morbidity and the need for general anesthesia^{1,2}. An alloplastic material avoids the need for a second surgical field but it should be safe, resorbable, able to maintain space, and cheap^{3,4}. Calcium sulfate (CaS) is highly biocompatible and it is one of the synthetic grafts with the longest clinical history (more than 100 years)⁵⁻¹². It has been utilized in treating periodontal disease, endodontic lesions, alveolar bone loss, and maxillary sinus augmentation^{3,4,13-18}. CaS has been used as a membrane to facilitate healing and to prevent the loss of other grafting materials¹⁹. When associated with other bone grafts it seems to have a favorable effect on osteogenesis^{2,20}. CaS rapidly resorbs and leaves a calcium phosphate lattice which promotes osteogenic activity^{21,22}. Ricci et al.²³ demonstrated that CaS induces new bone formation in dogs after 2 weeks and that it is almost completely resorbed after 1 month.

In previous studies we carried out a

genome wide screen of osteoblast-like cell line (MG-63) following treatment with CaS, using cDNA microarray. Several genes covering a broad range of functional activities, like signal transduction, differentiation, cell-cycle regulation, were significantly up-regulated²⁴. Then the genetic effect of CaS was studied in the same cell system at post-transcriptional level, with microRNA microarray²⁵. We identified miRNA that regulates the transduction of genes related to bone formation (TFIP11), skeletal development (MSX1, ADAMTS4, DLX5 FGFR1 COMP EN1 SHOX), cartilage remodeling (MATN1) and ossification (BMP1, BMP7, ALPL, PTH).

Because few reports analyze the genetic effects of CaS on stem cells²⁶, the expression of genes related to osteoblast differentiation were analyzed using cultures of mesenchymal stem cells derived from peripheral blood (PB-hMSCs) treated with CaS.

To investigate the osteogenic differentiation of PB-hMSCs, the quantitative expression of the mRNA of specific genes, like transcriptional factor (RUNX2), bone related genes (SPP1, COL1A1, COL3A1, BGLAP, ALPL, and FOSL1) and mesenchymal stem cells marker (ENG) were examined by means of real time Reverse Transcription-Polymerase Chain Reaction (real time RT-PCR).

MATERIALS AND METHODS

a) Stem preparation

PB-hMSCs were obtained for gradient centrifugation from peripheral blood of healthy anonymous volunteers, using the Acuspin System-Histopaque 1077 (Sigma Aldrich, Inc., St Louis, Mo, USA). Firstly, 30 ml of heparinized peripheral blood were added to the Acuspin System-Histopaque 1077 tube and centrifugated at 1000 x g for 10 minutes. After centrifugation the interface containing mononuclear cells was transferred in another tube, washed with PBS and centrifuged at 250 x g per 10 minutes. The enriched mononuclear pellets was resuspended in 10 ml of Alphamem medium (Sigma Aldrich, Inc., St Louis, Mo, USA) supplemented with antibiotics (Penicillin 100 U/ml and Streptomycin 100 micrograms/ml - Sigma, Chemical Co., St Louis, Mo, USA) and aminoacids (L-Glutamine - Sigma, Chemical Co., St Louis, Mo, USA). The cells were incubated at 37°C in a humidified atmosphere with 5% CO₂. Medium was changed after 24 hours. PB-hMSC were selected for adhesiveness and characterized as stem cells by immunofluorescence.

b) Immunofluorescence

Cells were washed with PBS for three times and fixed with cold methanol for 5 min at room temperature. After washing with PBS, cells were blocked with bovine albumin 3% (Sigma Aldrich, Inc., St Louis, Mo, USA) for 30 min at room temperature. The cells were incubated overnight sequentially at 4 °C with primary antibodies raised against CD105 1:200, mouse (BD Biosciences,

San Jose, CA, USA), CD73 1:200, mouse (Santa Cruz Biotechnology, Inc., Santa Cruz, CA, USA), CD90 1:200, mouse (Santa Cruz Biotechnology, Inc., Santa Cruz, CA, USA), CD34 1:200, mouse (Santa Cruz Biotechnology, Inc., Santa Cruz, CA, USA). They were washed with PBS and incubated for 1 h at room temperature with secondary antibody conjugated-Rodamine goat anti-mouse 1:200 (Santa Cruz Biotechnology, Inc., Santa Cruz, CA, USA). Subsequently, cells were mounted with the Vectashield Mounting Medium with DAPI (Vector Laboratories, Inc., Burlingame, CA, USA) and observed under a fluorescence microscope (Eclipse TE 2000-E, Nikon Instruments S.p.a., Florence, Italy).

c) Cell culture

PB-hMSCs at second passage were cultured in Alphamem medium (Sigma Aldrich, Inc., St Louis, Mo, USA) supplemented with 10% fetal calf serum, antibiotics (Penicillin 100 U/ml and Streptomycin 100 micrograms/ml - Sigma Aldrich, Inc., St Louis, Mo, USA) and aminoacids (L-Glutamine - Sigma Aldrich, Inc., St Louis, Mo, USA). The cells were incubated at 37°C in a humidified atmosphere with 5% CO₂.

For the assay, cells were collected and seeded at a density of 1x10⁵ cells/ml into 9 cm² (3ml) wells by using 0.1% trypsin, 0.02% EDTA in Ca⁺⁺ - and Mg - free Eagle's buffer for cell release.

One set of wells were added with CaS (Surgiplaster, Classimplant, Rome, Italy) at the concentration of 0.001 mg/ml. Another set of wells containing untreated cells were used as control. The medium was changed every 3 days.

After seven days, when cultures were sub-confluent, cells were processed for RNA extraction.

d) RNA processing

Reverse transcription to cDNA was performed directly from cultured cell lysate using the TaqMan Gene Expression Cells-to-Ct Kit (Ambion Inc., Austin, TX, USA), following manufacturer's instructions. Briefly, cultured cells were lysed with lysis buffer and RNA released in this solution. Cell lysate were reverse transcribed to cDNA using the RT Enzyme Mix and appropriate RT buffer (Ambion Inc., Austin, TX, USA). Finally the cDNA was amplified by real-time PCR using the included TaqMan Gene Expression Master Mix and the

specific assay designed for the investigated genes.

e) Real time PCR

Expression was quantified using real time RT-PCR. The gene expression levels were normalized to the expression of the housekeeping gene RPL13A and were expressed as fold changes relative to the expression of the untreated PB-hMSCs. Quantification was done with the delta/delta calculation method²⁷. Forward and reverse primers and probes for the selected genes were designed using primer express software (Applied Biosystems, Foster City, CA, USA) and are listed in **Table 1**.

All PCR reactions were performed in a 20 µl volume using the ABI PRISM 7500

(Applied Biosystems, Foster City, CA, USA). Each reaction contained 10 µl 2X TaqMan universal PCR master mix (Applied Biosystems, Foster City, CA, USA), 400 nM concentration of each primer and 200 nM of the probe, and cDNA. The amplification profile was initiated by 10-minute incubation at 95°C, followed by two-step amplification of 15 seconds at 95°C and 60 seconds at 60°C for 40 cycles. All experiments were performed including non-template controls to exclude reagents contamination. PCRs were performed with two biological replicates.

| Gene symbol | Gene name | Primer sequence (5'>3') | Probe sequence (5'>3') |
|-------------|-------------------------------------|---|------------------------------|
| SPP1 | osteopontin | F-GCCAGTTGCAGCCTTCTCA R-AAAAGCAAATCACTGCAATTCTCA | CCAAACGCCGACCAAGGAAAACCTCAC |
| COL1A1 | collagen type I alpha1 | F-TAGGGTCTAGACATGTTTCAGCTTTGT R-GTGATTGGTGGGATGTCTTCGT | CCTCTTAGCGGCCACCGCCCT |
| RUNX2 | runt-related transcription factor 2 | F-TCTACCACCCCGCTGTCTTC R-TGGCAGTGTCACTCATCTGAAATG | ACTGGGCTTCCTGCCATCACCGA |
| ALPL | alkaline phosphatase | F-CCGTGGCAACTCTATCTTTGG R-CAGGCCCATGCCATACAG | CCATGCTGAGTGACACAGACAAGAAGCC |
| COL3A1 | collagen, type III, alpha 1 | F-CCCACTATTATTTGGCACAACAG R-AACGGATCCTGAGTCACAGACA | ATGTTCCCATCTTGGTCAGTCCTATGCG |
| BGLAP | osteocalcin | F-CCCTCCTGCTTGGACACAAA R-CACACTCCTCGCCCTATTGG | CCTTTGCTGGACTCTGCACCGCTG |
| CD105 | endoglin | F-TCATCACCACAGCGGAAAAA R-GGTAGAGGCCAGCTGGAA | TGCACTGCCTCAACATGGACAGCCT |
| FOSL1 | FOS-like antigen 1 | F-CGCGAGCGGAACAAGCT R-GCAGCCCAGATTCTCATCTTC | ACTTCTGCAGGCGGAGACTGACAAAC |
| RPL13A | ribosomal protein L13 | F-AAAGCGGATGGTGGTTCCT R-GCCCCAGATAGGCAAACCTTC | CTGCCCTCAAGGTCGTGCGTCTG |

Table 1. Primer and probes used in real time PCR

RESULTS


PB-hMSCs were characterized by immunofluorescence. The cell surfaces were positive for mesenchymal stem cell markers, CD105, CD90 and CD73 and negative for marker of hematopoietic origin, CD34 (Figure 1). Transcriptional expressions of several osteoblast-related genes (RUNX2, SPP1, COL1A1, COL3A1, BGLAP, ALPL and FOSL1) and mesenchymal stem cells marker (ENG) were examined after 7 days of supplement treatment with CaS (0.001 mg/ml). CaS enhanced the expression of bone related genes like COL3A1, BGLAP and SPP1. The treatment did not affect the mRNA expression of ENG and FOSL1 that were similarly in both treated and untreated PB-hMSCs. RUNX2, COL1A1 and ALPL were decreased in the presence of CaS at day 7.

DISCUSSION

CaS is a highly biocompatible material⁵⁻¹². Solidified or crystallized CaS is very osteogenic *in vivo*. As the surface of CaS dissolves in body fluids, the calcium ions form calcium phosphate that re-precipitates on the surface forming an osteoblast “friendly” environment. How this “friendly” environment alters osteoblast activity to promote bone formation is incompletely understood. In order to get more inside how CaS acts on PB-hMSCs, changes in expression of bone related marker genes (RUNX2, SPP1, COL1A1, COL3A1, BGLAP, ALPL and FOSL1) and mesenchymal stem cells marker (ENG) were investigated by real-time RT-PCR.

Mesenchymal stem cells (MSCs) are defined as self-renewable, multipotent progenitors cells with the ability

Figure 1. PB-hMSCs by indirect immunofluorescence (Rodamine). Cultured cells were positive for the mesenchymal stem cell marker CD73 (b), CD90 (c), CD105 (d) and negative for the hematopoietic markers CD34 (a). Nucleuses were stained with DAPI. Original magnification x40.



to differentiate, under adequate stimuli, into several mesenchymal lineages, including osteoblasts²⁸. In our study, mesenchymal stem cells from human peripheral blood were isolated and characterized by morphology and immunophenotype. Isolated PB-hMSCs showed fibroblast-like morphology and were positive for MSCs surface molecules (CD90, CD105, CD73) and negative for markers of haematopoietic progenitors (CD34).

After 7 day of treatment with CaS the expression levels of osseodifferentiation genes were measured by relative quantification methods using real-time RT-PCR.

Two osteoblast-specific genes, SPP1 (an acid phosphoprotein involved in regulation of bone mineralization) and BGLAP (a bone specific protein in-

involved in mineralization and bone resorption), that are generally express by osteoblasts in the early stages of their differentiation²⁹ were significantly up-regulate in treated PB-hMSCs.

Expression of ENG and FOSL1 didn't have significant change in treated cells respect to control after 7 day of treatment with CaS. ENG (CD105) is a surface markers used to define a bone marrow stromal cell population capable of multilineage differentiation³⁰. This gene is a receptor for TGF- β 1 and - β 3³¹ and modulates TGF- β signaling by interacting with related molecules, such as TGF- β 1, - β 3, BMP-2, -7, and activin A. It is speculated that these members of the TGF- β superfamily are mediators of cell proliferation and differentiation and play regulatory roles in cartilage and bone formation³². The

disappearance of the CD105 antigen during osteogenesis suggests that this protein, like others in the TFG- β superfamily, is involved in the regulation of osteogenesis³³. FOSL1 that encodes for Fra-1, a component of the dimeric transcription factor activator protein-1 (AP-1), which is composed mainly of Fos (c-Fos, FosB, Fra-1 and Fra-2) and Jun proteins (c-Jun, JunB and JunD). AP-1 sites are present in the promoters of many developmentally regulated osteoblast genes, including alkaline phosphatase and collagen I. McCabe et al.³⁴ demonstrated that differential expression of Fos and Jun family members could play a role in the developmental regulation of bone-specific gene expression and, as a result, may be functionally significant for osteoblast differentiation. Kim et al.³⁵ studying the effect of a new anabolic agents that stimulate bone formation, find that this gene is activated in the late stage of differentiation, during the calcium deposition.

In our study FOSL1 and ENG were weakly down- and up-regulated respectively, probably because cells were at early stage of differentiation. RUNX2 and ALPL were down regulated in treated cells respect to control after 7 day of treatment with CaS. RUNX2 is the most




Figure 2. Gene expression analysis of PB-hMSCs after 7 days of treatment with CaS.

specific osteoblast transcription factor and is a prerequisite for osteoblast differentiation and consequently mineralization. This result is comparable with data reported by Kim et al.³⁶ They showed that there was no BMP-2-mediated up-regulation of RUNX2 mRNA expression at days 3 or 7³⁶.

Alkaline phosphatase regulates mineralization of bone matrix. Several studies demonstrated that the potency of individual substances to induce alkaline phosphatase varies in a species-dependent manner. Glucocorticoids such as dexamethasone are potent inducers in human and rat stromal cells, but they have no effect on alkaline phosphatase activity in mouse stromal cells^{37,38}. On the contrary, bone morphogenetic proteins (BMPs) are potent inducers of osteogenesis in both mouse and rat bone marrow stromal cells³⁹ but Diefenderfer et al showed that BMP-2 alone is a poor osteoblast inducer in human marrow derived stromal cells⁴⁰.

CaS also modulates the expression of genes encoding for collagenic extracellular matrix proteins like collagen type 1 α 1 (COL1A1). Collagen type1 is the most abundant in the human organism⁴¹. In our study COL1A1 is significantly down expressed as compared to the control when exposed to CaS probably because this gene is activated in the late stage of differentiation and are related to extracellular matrix synthesis.

On the contrary, COL3A1 was up-regulated in treated cells respect to control. COL3A1 encodes the pro-alpha1 chains of type III collagen, a fibrillar collagen that is found in extensible connective tissues⁴².

CONCLUSIONS

The present study shows the effect of CaS on PB-hMSCs in the early differentiation stages: CaS is an inducer of osteogenesis on human stem cells as demonstrated by the activation of bone related genes: osteopontin (SPP1), osteocalcin (BGLAP) and COL3A1. Moreover, we chose to perform the experiment after 7 days in order to get information on the early stages of stimulation. It is our understanding, therefore, that more investigations with different time points are needed in order to get a global comprehension of the molecular events related to CaS action. The reported model is useful to investigate the effects of different substances on stem cells.

Acknowledgments

This work was supported by FAR from the University of Ferrara (FC), Ferrara, Italy, and from Regione Emilia Romagna, Programma di Ricerca Regione Università, 2007–2009, Area 1B: Patologia osteoarticolare: ricerca pre-clinica e applicazioni cliniche della medicina rigenerativa, Unità Operativa n. 14.

REFERENCES

1. Guarnieri R, Bovi M. Maxillary sinus augmentation using prehardened calcium sulfate: a case report. *Int J Periodontics Restorative Dent* 2002;22:503-8.
2. Macneill SR, Cobb C M, Rapley J W, Glaros G, Spencer P. In vivo comparison of synthetic osseous graft materials. A preliminary study. *J Clin Periodontol* 1999;26:239-45.
3. Pecora G, De Leonardis D, Ibrahim N, Bovi M, Cornelini R. The use of calcium sulphate in the surgical treatment of a 'through and through' periradicular lesion. *Int Endod J* 2001;34:189-97.
4. Yoshikawa G, Murashima Y, Wadachi R, Sawada N, Suda H. Guided bone regeneration

- (GBR) using membranes and calcium sulphate after apicectomy: a comparative histomorphometrical study. *Int Endod J* 2002;35:255-63.
5. Alderman N. Sterile plaster of Paris as an implant in the infrabony environment: a preliminary study. *J Periodontol* 1969;40:11-13.
 6. Alexander DJ, Manson NA, Mitchell MJ. Efficacy of calcium sulfate plus decompression bone in lumbar and lumbosacral spinal fusion: preliminary results in 40 patients. *Can J Surg* 2001;44:262-6.
 7. Gitelis S, Piasecki P, Turner T, Haggard W, Charters J, Urban R. Use of a calcium sulfate-based bone graft substitute for benign bone lesions. *Orthopedics* 2001;24:162-6.
 8. Peltier LF. The use of plaster of Paris to fill large defects in bone. *Am J Surg* 1959;97:311.
 9. Pietrzak WS, Ronk R. Calcium sulfate bone void filler: a review and a look ahead. *J Craniofac Surg* 2000;11:327-33; discussion 34.
 10. Shaffer CD, App GR. The use of plaster of paris in treating infrabony periodontal defects in humans. *J Periodontol* 1971;42:685-90.
 11. Sottosanti JS. Calcium sulfate-aided bone regeneration: a case report. *Periodontol Clin Invest* 1995;17:10-5.
 12. Coetzee a S. Regeneration of bone in the presence of calcium sulfate. *Arch Otolaryngol* 1980;106:405-9.
 13. Bier SJ, Sinensky MC. The versatility of calcium sulfate: resolving periodontal challenges. *Compend Contin Educ Dent* 1999;20:655-61; quiz 62.
 14. Kim CK, Kim HY, Chai JK, Cho KS, Moon IS, Choi SH, Sottosanti JS, Wikesjo U M. Effect of a calcium sulfate implant with calcium sulfate barrier on periodontal healing in 3-wall intrabony defects in dogs. *J Periodontol* 1998;69:982-8.
 15. Orsini M, Orsini G, Benlloch D, Aranda JJ, Lazaro P, Sanz M, De Luca M, Piattelli A. Comparison of calcium sulfate and autogenous bone graft to bioabsorbable membranes plus autogenous bone graft in the treatment of intrabony periodontal defects: a split-mouth study. *J Periodontol* 2001;72:296-302.
 16. Pecora G, Baek SH, Rethnam S, Kim S. Barrier membrane techniques in endodontic microsurgery. *Dent Clin North Am* 1997;41:5817.
 17. Pecora GE, De Leonardis D, Della Rocca C, Cornelini R, Cortesini C. Short-term healing following the use of calcium sulfate as a grafting material for sinus augmentation: a clinical report. *Int J Oral Maxillofac Implants* 1998;13:866-73.
 18. Kelly CM, Wilkins RM, Gitelis S, Hartjen C, Watson JT, Kim P T. The use of a surgical grade calcium sulfate as a bone graft substitute: results of a multicenter trial. *Clin Orthop Relat Res* 2001;42-50.
 19. Maragos P, Bissada NF, Wang R, Cole BP. Comparison of three methods using calcium sulfate as a graft/barrier material for the treatment of Class II mandibular molar furcation defects. *Int J Periodontics Restorative Dent* 2002;22:493-501.
 20. Al Ruhaimi KA. Effect of adding resorbable calcium sulfate to grafting materials on early bone regeneration in osseous defects in rabbits. *Int J Oral Maxillofac Implants* 2000;15:859-64.
 21. Damien CJ, Ricci JL, Christel P, Alexander H, Patat JL. Formation of a calcium phosphate-rich layer on absorbable calcium carbonate bone graft substitutes. *Calcif Tissue Int* 1994;55:151-8.
 22. Tay BK, Patel VV, Bradford D S. Calcium sulfate- and calcium phosphate-based bone substitutes. Mimicry of the mineral phase of bone. *Orthop Clin North Am* 1999;30:615-23.
 23. Ricci JL, Nadkarni P, Hawkins M, Turner J, Rosenblum S, Brezenoff L, De Leonardis D, Pecora G. Biological mechanisms of calcium sulfate replacement by bone ed. J.E. Davies. 2000, Toronto: Em Squared Incorporated. 332-44.
 24. Carinci F, Piattelli A, Stabellini G, Palmieri A, Scapoli L, Laino G, Caputi S, Pezzetti F. Calcium sulfate: analysis of MG63 osteoblast-like cell response by means of a microarray technology. *J Biomed Mater Res B Appl Biomater* 2004;71:260-7.
 25. Palmieri A, Pezzetti F, Brunelli G, Scapoli L, Lo Muzio L, Scarano A, Martinelli M, Carinci F. Calcium sulfate acts on the miRNA of MG63E osteoblast-like cells. *J Biomed Mater Res B Appl Biomater* 2008;84:369-74.
 26. Lazary A, Balla B, Kosa JP, Bacsi K, Nagy Z, Takacs I, Varga PP, Speer G, Lakatos P. Effect of gypsum on proliferation and differentiation of MC3T3-E1 mouse osteoblastic cells. *Biomaterials* 2007;28:393-9.
 27. Livak KJ, Schmittgen TD. Analysis of relative gene expression data using real-time quantitative PCR and the 2(-Delta Delta C(T)) Method. *Methods* 2001;25:402-8.
 28. Alhadlaq A, Mao JJ. Mesenchymal stem cells: isolation and therapeutics. *Stem Cells Dev* 2004;13:436-48.
 29. Lian Jb a J. Bone formation: maturation and functional activities of osteoblast lineage cells, in *Primer on the metabolic bone diseases and disorders on mineral metabolism*, M.J. Favus, Editor. 2003, The American Society for Bone and Mineral Research: Washington p. 13-28.
 30. Pittenger MF, Mackay aM, Beck SC, Jaiswal RK, Douglas R, Mosca JD, Moorman MA, Simonetti DW, Craig S, Marshak DR. Multilineage potential of adult human mesenchymal stem cells. *Science* 1999;284:143-7.
 31. Barry FP, Boynton RE, Haynesworth S, Murphy JM, Zaia J. The monoclonal antibody SH-2, raised against human mesenchymal stem cells, recognizes an epitope on endoglin (CD105). *Biochem Biophys Res Commun* 1999;265:134-9.
 32. Jakob M, Demarteau O, Schafer D, Hintermann B, Dick W, Heberer M, Martin I. Specific growth factors during the expansion and redifferentiation of adult human articular chondrocytes enhance chondrogenesis and cartilaginous tissue formation in vitro. *J Cell Biochem* 2001;81:368-77.
 33. Haynesworth SE, Baber MA, Caplan a I. Cell surface antigens on human marrow-derived mesenchymal cells are detected by monoclonal antibodies. *Bone* 1992;13:69-80.
 34. McCabe LR, Banerjee C, Kundu R, Harrison RJ, Dobner PR, Stein JL, Lian JB, Stein GS. Developmental expression and activities of specific fos and jun proteins are functionally related to osteoblast maturation: role of Fra-2 and Jun D during differentiation. *Endocrinology* 1996;137:4398-408.
 35. Kim JM, Lee SU, Kim YS, Min YK, Kim SH. Bone morphogenetic protein-2 stimulates osteoblast differentiation via coordinating activation of MAP kinases and transcription factors. *J Cell Biochem* 2008;104:1906-17.
 36. Kim IS, Song YM, Cho TH, Park YD, Lee KB, Noh I, Weber F, Hwang SJ. In vitro response of primary human bone marrow stromal cells to recombinant human bone morphogenetic protein-2 in the early and late stages of osteoblast differentiation. *Dev Growth Differ* 2008;
 37. Leboy PS, Beresford JN, Devlin C, Owen ME. Dexamethasone induction of osteoblast mRNAs in rat marrow stromal cell cultures. *J Cell Physiol* 1991;146:370-8.
 38. Beresford JN, Joyner CJ, Devlin C, Triffitt JT. The effects of dexamethasone and 1,25-dihydroxyvitamin D3 on osteogenic differentiation of human marrow stromal cells in vitro. *Arch Oral Biol* 1994;39:941-7.
 39. Balk ML, Bray J, Day C, Epperly M, Greenberger J, Evans CH, Niyibizi C. Effect of rhBMP-2 on the osteogenic potential of bone marrow stromal cells from an osteogenesis imperfecta mouse (oim). *Bone* 1997;21:7-15.
 40. Diefenderfer DL, Osyczka aM, Garino JP, Leboy PS. Regulation of BMP-induced transcription in cultured human bone marrow stromal cells. *J Bone Joint Surg Am* 2003;85-A Suppl 3:19-28.
 41. Suuriniemi M, Kovanen V, Mahonen A, Alen M, Wang Q, Lyytikainen A, Cheng S. COL1A1 Sp1 polymorphism associates with bone density in early puberty. *Bone* 2006;39:591-7.
 42. Chan TF, Poon A, Basu A, Addleman NR, Chen J, Phong A, Byers P H, Klein TE, Kwok PY. Natural variation in four human collagen genes across an ethnically diverse population. *Genomics* 2008;91:307-14.

World Leader

Safescraper[®]

Cortical Bone Collector **TWIST**



Versione dritta

Versione curve

Garantisci ai tuoi pazienti il Gold Standard

L'osso autologo corticale è sicuro, osteogenico e può essere utilizzato in combinazione con altri sostituti ossei. Lo strumento è monouso, blisterato sterile singolarmente

Il migliore strumento di prelievo osseo: 45-72% di vitalità cellulare

La tecnica di prelievo manuale preserva la componente cellulare dell'innesto. I riccioli corticali ottenuti con la lama del Safescraper[®] presentano cellule ossee vitali e ben preservate

Leader in qualità

La lama ottenuta attraverso tecnologie avanzatissime è il cuore dello strumento e rende Safescraper[®] unico nel mercato

Zaffe D, D'Avenia F. A novel bone scraper for intraoral harvesting: a device for filling small bone defects. Clin Oral Implants Res 2007; 18 (4): 525-33

Guida L, Annunziata M, Belardo S, Farina R, Scabbia A, Trombelli L. Effect of autogenous cortical bone particulate in conjunction with enamel matrix derivative in the treatment of periodontal intraosseous defects. J Periodontol 2007; 78 (2): 231-8

Trombelli L, Farina R, Marzola A, Ibro A, Calura G. GBR and autogenous cortical bone particulate by bone scraper for alveolar ridge augmentation: a 2-case report. Int J Oral Maxillofac Implants 2008 Jan-Feb; 23 (1): 111-6



CE
0123



C.G.M. S.p.A.
DIVISIONE MEDICALE META
Via E. Villa, 7 - 42100 Reggio Emilia
Tel.: 0522.50.23.11 - Fax: 0522.50.23.33 - www.metahosp.com

Zimmer...non solo Impianti !

Puros[®]
Blocco per
Alloinnesto



Membrana Zimmer
Socket Repair[®]

Puros[®]
Particolato di Osso
Corticale/Spongioso per Alloinnesto



CopiOs,
Membrana Pericardica



Tapered
Screw-Vent



Spline



Swiss Plus



One Piece



zimmer
dental

Collaborare con la natura per riempire un vuoto!

Zimmer Dental Italy srl Via Matteotti, 98 - 31029 Vittorio Veneto (TV) - Italy
Tel.0438 555573 r.a. Fax 0438 553181 zimmerdental.italy@zimmer.com